

Test Strip Detection

Please refer to information provided on a specific cup, details here are general, but products may vary.

Common Name	Strip	Type	Min Detection	Max Detection	Cut-off Levels
Amphetamines	AMP AMP	AMP	2-7 hours	2-4 Days	1000 ng/mL
Barbiturates	BAR BAR	BAR	2-4 hours	2-10 Days	300 ng/mL
Buprenorphine	BUP BUP	BUP	2-4 hours	2-6 Days	10 ng/mL
Benzodiazepines	BZO BZO	BZO	2-7 hours	1-6 Weeks	300 ng/mL
Cocaine	COC COC	COC	1-4 hours	2-7 Days	300 ng/mL
Methadone	EDDP-2 EDDP-2	EDDP-2	1 day	30 days	300 ng/mL
Alcohol	ETG ETG	EtG	1-3 hours	4-5 Days	500 ng/mL
Fentanyl	FEN FEN	FYL/FEN	1 hour	3 Days	200 ng/mL
Gabapentin	GAB GAB	GABA	4-6 hours	-	2000 ng/mL
Spice/K2	K2 K2	K2	1-2 hours	1-2 Days	50 ng/mL
Kratom	KRA KRA	KRA	1-5 days	24 hours	250 ng/mL
Ecstasy	MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA MDMA	MDMA	2-7 hours	2-4 Days	500 ng/mL
Methamphetamine	MET MET	MET/mAMP	2-7 hours	2-5 Days	1000 ng/mL
Methadone	MTD MTD	MTD	3-8 hours	1-3 Days	300 ng/mL
Opiates	OPI OPI	OPI	2 hours	2-4 Days	2000 ng/mL
Opiates 300	MOP MOP	OPI 300	2 hours	2-4 Days	300 ng/mL
Oxycodone	OXY OXY	OXY	1-3 hours	1-2 Days	100 ng/mL
Phencyclidine	PCP PCP	PCP	4-6 hours	7-14 Days	25 ng/mL
Propoxyphene	PPX PPX	PPX	1-6 hours	7-10 Days	300 ng/mL
Tricyclic Antidepressants	TCA TCA	TCA	8-12 hours	2-7 Days	1000 ng/mL
Marijuana	THC THC	THC	2 hours	Up to 14 Days	50 ng/mL
Tramadol	TRA TRA	TML/TRA	1-5 days	-	100 ng/mL